



Otolift Modul-Air

A revolutionary stairlift



Content

Revolutionary stairlift	4
The advantages	8
Rail cross section	12
Reverse Drive Technology	13
Flexible comfort	14
Intuitive and easy to use	16
Pick your upholstery and paint colour	20
Arrange an assessment	22



Otolift Modul-Air

Revolutionary stairlift

Otolift Stairlifts released a new stairlift onto the market this year, the Otolift Modul-Air. The stairlift with the thinnest single rail in the world is now available for modular delivery.

In 2011 we introduced the Otolift One, a stairlift on a single rail. This stairlift received the Red Dot Design Award. In 2014 the Otolift Air was created as successor to the Otolift One. The same design, but a 25 percent thinner rail. Now, two years later, it is time for a new revolution: the Otolift Modul-Air. The stairlift with the thinnest single rail in the world and a beautiful design offers many advantages to you as a customer.



reddot design award



Otolift Modul-Air

The revolution
amongst the stairlifts

Due to modular technology it is possible to keep all rail parts in stock. This can not only save costs, but also reduce delivery time.



Otolift Modul-Air

The advantages

What are the advantages of the revolutionary Otolift Modul-Air?

The Modul-Air has the thinnest single rail in the world, consequently it blends in with your interior. The stairlift is always provided with a powered folding footrest and can be provided with an automated swivel seat.



Otolift Modul-Air

The 6 advantages

1. No customisation, better value
2. Thinnest single rail in the world
3. Super compact seat
4. Automated swivel seat
5. Powered folding footrest
6. Fixed to the stairs

Unique!
Not costum,
so inexpensive



Otolift Modul-Air Rail cross-section

This makes it particularly clear: the predecessor (Otolift One) and the new rail system of the Otolift Modul-Air next to each other.



Otolift One (predecessor)
ø 8 cm in real size



Otolift Modul-Air
ø 6 cm in real size

Otolift Modul-Air Reverse Drive Technology

Thanks to this smart technology the lift can travel partly or even entirely backwards. This prevents you from hitting the wall or handrail with your knees.



Otolift Modul-Air

Flexible convenience

The Otolift Modul-Air is very suitable for installation alongside the broad or narrow side of your stairs. The thin rail of the Otolift Modul-Air is also hardly noticeable, thus your stairs hardly look any different.

Broad, narrow or circular stairs, the Otolift Modul-Air easily fits the most demanding stairs and easily moves with all turns. Furthermore, the rail takes up so little space that the stairs can also be easily used by foot. In this way, it is guaranteed that not only you, but also your family members or guests, can always safely use the stairs.



Automatically folding footrest

Equipped with an automatically folding footrest as a standard. This way you do not need to bend down to fold the footrest.

Automated swivel seat

Equipped with an automated swivel seat, for a smooth ride and getting on and off safely. (Optional)

Reverse Drive Technology

Travelling backwards is also an option, providing more seating comfort. Also this prevents you from hitting walls or hand rails with your knees. (Optional)



Otolift Modul-Air

Intuitive and easy to use

Your Otolift was designed for convenience. The use is very intuitive and will soon become second nature.



Step 1

You move the chair towards you using the wireless controller.



Step 2

You unfold the chair lightly and easily.



Step 3

You sit down on the comfortable seat.



Step 4

Operate the control and the chair starts moving quietly and smoothly.



Step 5

The chair quietly continues its way upstairs and stops automatically.



Step 6

The seat turns automatically (depending on model) after which you step off safely.

Otolift Modul-Air

Intuitive and easy to use



Safe

The safety belt is easy to use and can be operated with one hand.



Smart

The stairlift keeps functioning in case of power failure.



Intelligent

Your Otolift will stop when detecting an obstacle. Automatic, soft and safe.



User-friendly

The control lever is robust and stable enough to be used as a support.



Space saving

The Otolift can be folded quickly. This saves space and leaves the stairs accessible.



Short start

Because of the short start, the passage or door at the bottom of the stairs remains clear.

Otolift Modul-Air

Pick your upholstery and paint colour

The new Otolift Modul-Air does not only persuade by its exclusive design, but also because of the options available, that allow an element of personalisation to your lift. The Otolift is available with a seat in your favourite colour for a perfect match with your interior. Also in terms of materials there are various options. You choose and we make it happen.

Upholstery seat and backrest

Fabric



Grey



Cream



Red



Brown

Synthetic leather



Hazelnut



Pearl white




Blue




Black

Leather



Taupe




Chestnut

Paint work



Creme white
RAL 9001



Brown
RAL 8025



Arrange an Assessment for an Otolift

We have already provided a lot of information about our unique stairlifts. But what will be the cost of your Otolift?

Unfortunately as every staircase is different, and you have your own preferences and wishes, we cannot give a simple answer to this question without assessing your staircase and individual needs in person.

Free and non-binding staircase evaluation.

Our experienced advisors work all across the country and will be happy to visit you to give personal advice. Of course this is done free of charge and entirely voluntarily. Our adviser will evaluate your staircase and listen well to your individual wishes and preferences. He will give you transparent advice about which Otolift best fits you and along which side of the stairs it could best be installed.

At that moment, for us it is possible to determine the actual price of your Otolift.

Of course we hope that you will consider us when taking your next step: a non-binding staircase assessment!





UPSTAIRS

Upstairs Sdn Bhd (918460-P)

C-5-1, Pusat Perniagaan Bukit Serdang,
Jalan BS 14/3 Taman Bukit Serdang,
Seksyen 14, 43300 Seri Kembangan Selangor
Malaysia

Contact No.: +603 - 9548 0100
Email: sales@upstairs.com.my
Website: www.upstairs.com.my

